

AccuVein

AccuVein[®] AV300 Vein illumination device

- . Easy to use and transport
- . Builds confidence of phlebotomists
- . Increases staff productivity and decreases costs
- . Designed to be non-contact, no sterilization between uses

The AccuVein AV300 vein illumination device allows phlebotomists and nurses to quickly and accurately locate donors and patients veins beneath the skin.

In blood centers

A successful venipuncture is an essential element in ensuring donor satisfaction and a successful blood draw. Deferrals and incomplete blood draws are inefficient and costly to operations.

In hospitals

Whether drawing blood specimens or establishing IV access, venipuncture is the most commonly performed invasive medical procedure. It is also one that frequently requires additional time and multiple attempts to be successful.

Find out today how AccuVein can help keep operations flowing.

AccuVein[®] AV300

Vein illumination device

The AccuVein AV300 works using infrared light technology which is absorbed by hemoglobin in the blood. When the AccuVein is held about 18 cm above the skin, veins appear noticeably different than the surrounding tissue. The vasculature shows up clearly on the skin's surface, aiding in vein location for accurate venipuncture.

Easy to Use and Transport

The AccuVein AV300 offers a simple point and click design that requires no pre-use calibration or adjustments. The device fits in your hand, weighs only 283 gr and comes with a rechargeable battery that allows for easy use in mobile or fixed site environments. Also the AccuVein AV300 has been designed for non-contact operation.

Benefits of the AccuVein AV300 in blood centers

- . Helps visualize the best vein for each procedure and maximizes efficiency and revenue by minimizing loss of component revenue due to vasculature deferrals and QNS.
- . Enhances perception of a positive donation experience for donors by improving vein visualization and donor assessment for a successful first stick (mitigates fear of needles).

Benefits of the AccuVein AV300 in hospitals

. Increases staff productivity: a missed vein requires the procedure to be repeated; productivity can be further impacted should the phlebotomist have to escalate the procedure to specialized staff.

- Provides more timely treatment/diagnosis: blood tests that facilitate diagnosis and patient management may be delayed because of difficult venipuncture.
- . Decreases patient discomfort: patients may endure unnecessary needle sticks.
- . Decreases stress for both patients and staff.
- . Helps to achieve specific cost reduction in: unnecessary PICC lines, unnecessary use of IV start kits.
- . Decreases the risk of infections (catheter related blood infections).



For more information, please contact your local Sales Representative.

Popular accessories for the AccuVein AV300



Belt carry case

Powered hands-free kit



Unpowered hands-free kit



Fenwal Europe 2 - 4 Rue Edouard Belin

B - 1435 Mont Saint Guibert Belgium www.fenwalinc.com Fenwal is a registered trademark of Fenwal Inc. in the United States and other countries. ACCUVEIN is a trademark of AccuVein, LLC.